



Maxim Biotech, Inc.

Certificate of Analysis

Cat. No.: IH-60025 (MMP2-6001-B)

Biotinylated Long DNA Probe for Human MMP2

 8 µg / tube

Form: Liquid 8 µg biotin-labelled DNA in 100 µl ddH₂O

Tm: > 90 °C (Calculated according to the nearest neighborhood method)

Size: 191bp. This biotin-conjugated probe was produced by PCR primers (Maxim's product Catalog # SP-10514) under PCR* conditions using Human cDNA and biotin-dUTPs & biotin-dATPs.

Purity: The biotin-conjugated probe was purified by cartridge and showed a single band on gel.

Storage and Stability: Recommend to store at -20 °C, stable for one year. It is stable for one month when store at 4 °C.

Use: The product is able to detect a fragment of **Human MMP2** in in-situ hybridization when used together with Maxim's in-situ hybridization /detection kit (IHD-0050 & IHD0052) and/or in Northern hybridization.

Sequences: (Alignment on database: Genbank, **BC002576**)

References: Strausberg,R. Submitted (05-FEB-2001) National Institutes of Health, Mammalian Gene Collection (MGC), Cancer Genomics Office, National Cancer Institute, 31 Center Drive, Room 11A03, Bethesda, MD 20892-2590, USA

Lot. No.: _____

Date of Production: _____

Operator: _____

PCR*: The Polymerase Chain Reaction (PCR) process is covered by patents owned by Hoffmann-La Roche. Use of the PCR Process requires a license. This product does not convey a license to use the PCR process.

Maxim Biotech Inc.

780 Dubuque Avenue, South San Francisco, CA 94080 USA
Tel: 650-871-1919 Fax: 650-871-2857 E-mail: mbi@maximbio.com