



Maxim Biotech, Inc.

Certificate of Analysis

Cat. No.: IH-60048 (FAS-6010-B)
Biotinylated Long DNA Probe for Human FasL Gene 8 µg / tube

Form: Liquid 8 µg biotin-labelled DNA in 100 µl ddH₂O

Tm: > 90 °C (Calculated according to the nearest neighborhood method)

Size: 250bp. This biotin-conjugated probe was produced by PCR primers (Maxim's product Catalog # SP-10215) under PCR* conditions using Human cDNA and biotin-dUTPs & biotin-dATPs.

Purity: The biotin-conjugated probe was purified by cartridge and showed a single band on gel.

Storage and Stability: Recommend to store at -20 °C, stable for one year. It is stable for one month when store at 4 °C.

Use: The product is able to detect a fragment of **Human FasL Gene** in in-situ hybridization when used together with Maxim's in-situ hybridization /detection kit (IHD-0050 & IHD0052) and/or in Northern hybridization.

Sequences: (Alignment on database: Genbank, **D38122, U08137**)

References: Mita,E., et al. Biochem. Biophys. Res. Commun. 204 (2), 468-474 (1994)

Lot. No.: _____

Date of Production: _____

Operator: _____

PCR*: The Polymerase Chain Reaction (PCR) process is covered by patents owned by Hoffmann-La Roche. Use of the PCR Process requires a license. This product does not convey a license to use the PCR process.

Maxim Biotech Inc.

780 Dubuque Avenue, South San Francisco, CA 94080 USA
Tel: 650-871-1919 Fax: 650-871-2857 E-mail: mbi@maximbio.com