



Maxim Biotech, Inc.

Certificate of Analysis

Cat. No.: IH-60050 (GMC-6001-B)
Biotinylated Long DNA Probe for Human GM-CSF 8 µg / tube

Form: Liquid 8 µg biotin-labelled DNA in 100 µl ddH₂O

Tm: > 90 °C (Calculated according to the nearest neighborhood method)

Size: 424bp. This biotin-conjugated probe was produced by PCR primers (Maxim's product Catalog # SP-10249) under PCR* conditions using Human cDNA and biotin-dUTPs & biotin-dATPs.

Purity: The biotin-conjugated probe was purified by cartridge and showed a single band on gel.

Storage and Stability: Recommend to store at -20 °C, stable for one year. It is stable for one month when store at 4 °C.

Use: The product is able to detect a fragment of **Human GM-CSF** in in-situ hybridization when used together with Maxim's in-situ hybridization /detection kit (IHD-0050 & IHD0052) and/or in Northern hybridization.

Sequences: (Alignment on database: Genbank, **NM_000758**)

References: Wong,G.G., et al. Science 228 (4701), 810-815 (1985)

Lot. No.:

Operator:

PCR*: The Polymerase Chain Reaction (PCR) process is covered by patents owned by Hoffmann-La Roche. Use of the PCR Process requires a license. This product does not convey a license to use the PCR process.

Maxim Biotech Inc.

780 Dubuque Avenue, South San Francisco, CA 94080 USA
Tel: 650-871-1919 Fax: 650-871-2857 E-mail: mbi@maximbio.com