



# Maxim Biotech, Inc.

## Certificate of Analysis

**Cat. No.:** IH-60056 (HPV-6007)

**Biotinylated Long DNA Probe for HPV Type 6**        8   µg / tube

**Form:** Liquid 8 µg biotin-labelled DNA in 100 µl ddH<sub>2</sub>O

**Tm:** > 90 °C (Calculated according to the nearest neighborhood method)

**Size:** 263bp. This biotin-conjugated probe was produced by PCR primers (Maxim's product Catalog # SP-10298) under PCR\* conditions using Human cDNA and biotin-dUTPs & biotin-dATPs.

**Purity:** The biotin-conjugated probe was purified by cartridge and showed a single band on gel.

**Storage and Stability:** Recommend to store at -20 °C, stable for one year. It is stable for one month when store at 4 °C.

**Use:** The product is able to detect a fragment of **HPV Type 6** in in-situ hybridization when used together with Maxim's in-situ hybridization /detection kit (IHD-0050 & IHD0052) and/or in Northern hybridization.

**Sequences:** (Alignment on database: Genbank, **X00203**)

**References:** Schwarz,E., et al. EMBO J. 2 (12), 2341-2348 (1983)

**Lot. No.:** \_\_\_\_\_

**Date of Production:** \_\_\_\_\_

**Operator:** \_\_\_\_\_

PCR\*: The Polymerase Chain Reaction (PCR) process is covered by patents owned by Hoffmann-La Roche. Use of the PCR Process requires a license. This product does not convey a license to use the PCR process.

### Maxim Biotech Inc.

780 Dubuque Avenue, South San Francisco, CA 94080 USA  
Tel: 650-871-1919 Fax: 650-871-2857 E-mail: mbi@maximbio.com