



Maxim Biotech, Inc.

Certificate of Analysis

Cat. No.: IH-60061 (HPV-6001-B)
Biotinylated Long DNA Probe for *H. pylori*/Cag A 8 µg / tube

Form: Liquid 8 µg biotin-labelled DNA in 100 µl ddH₂O

Tm: > 90 °C (Calculated according to the nearest neighborhood method)

Size: 349bp. This biotin-conjugated probe was produced by PCR primers (Maxim's product Catalog # SP-10304) under PCR* conditions using Human cDNA and biotin-dUTPs & biotin-dATPs.

Purity: The biotin-conjugated probe was purified by cartridge and showed a single band on gel.

Storage and Stability: Recommend to store at -20 °C, stable for one year. It is stable for one month when store at 4 °C.

Use: The product is able to detect a fragment of ***H. pylori*/Cag A** in in-situ hybridization when used together with Maxim's in-situ hybridization /detection kit (IHD-0050 & IHD0052) and/or in Northern hybridization.

Sequences: (Alignment on database: Genbank, **L11714**)

References: Tummuru, M.K., et al. Infect. Immun. 61 (5), 1799-1809 (1993)

Lot. No.: _____

Date of Production: _____

Operator: _____

PCR*: The Polymerase Chain Reaction (PCR) process is covered by patents owned by Hoffmann-La Roche. Use of the PCR Process requires a license. This product does not convey a license to use the PCR process.

Maxim Biotech Inc.

780 Dubuque Avenue, South San Francisco, CA 94080 USA
Tel: 650-871-1919 Fax: 650-871-2857 E-mail: mbi@maximbio.com