



Maxim Biotech, Inc.

Certificate of Analysis

Cat. No.: IH-60063 (HSV-6007)

Biotinylated Long DNA Probe for HSV I/II DNA Pol. 8 µg / tube

Form: Liquid 8 µg biotin-labelled DNA in 100 µl ddH₂O

Tm: > 90 °C (Calculated according to the nearest neighborhood method)

Size: 92bp. This biotin-conjugated probe was produced by PCR primers (Maxim's product Catalog # HSV-1012/1013) under PCR* conditions using HSV Type II DNA and biotin-dUTPs & biotin-dATPs.

Purity: The biotin-conjugated probe was purified by cartridge and showed a single band on gel.

Storage and Stability: Recommend to store at -20 °C, stable for one year. It is stable for one month when store at 4 °C.

Use: The product is able to detect a fragment of **HSV I/II DNA Pol.** in in-situ hybridization when used together with Maxim's in-situ hybridization /detection kit (IHD-0050 & IHD0052) and/or in Southern hybridization.

Sequences: (Alignment on database: Genbank, [NC_001806](#), [NC_001798](#))

References: McGeoch,D.J., et al. J. Mol. Biol. 181 (1), 1-13 (1985)

Lot. No.: _____

Date of Production: _____

Operator: _____

PCR*: The Polymerase Chain Reaction (PCR) process is covered by patents owned by Hoffmann-La Roche. Use of the PCR Process requires a license. This product does not convey a license to use the PCR process.

Maxim Biotech Inc.

780 Dubuque Avenue, South San Francisco, CA 94080 USA
Tel: 650-871-1919 Fax: 650-871-2857 E-mail: mbi@maximbio.com