

Rapid Test Form for Swift HIV Recent Infection Assay (RIA)

Test Name:				Sample Type (circle):	Finger Prick	Venous Blood	Plasma	Serum	
Kit Lot#		Kit Expiry Date		Site/ Facility Name			Test Name		
No.	Test Date	Sample ID	Visual Result Mark "√" if line is present			Recency Interpretation All 3 lines = LT; C & T lines = Recent; Only C line = Neg; Others = Invalid (Please circle one)			
			Control line	Test line	Long Term (LT) line				
		QC - Long Term	√	√	√	(LT)	Recent	Neg	Invalid
		QC - Recent	√	√		LT	(Recent)	Neg	Invalid
		QC - Negative	√			LT	Recent	(Neg)	Invalid
1.						LT	Recent	Neg	Invalid
2.						LT	Recent	Neg	Invalid
3.						LT	Recent	Neg	Invalid
4.						LT	Recent	Neg	Invalid
5.						LT	Recent	Neg	Invalid
6.						LT	Recent	Neg	Invalid
7.						LT	Recent	Neg	Invalid
8.						LT	Recent	Neg	Invalid
9.						LT	Recent	Neg	Invalid
10.						LT	Recent	Neg	Invalid
11.						LT	Recent	Neg	Invalid
12.						LT	Recent	Neg	Invalid
13.						LT	Recent	Neg	Invalid
14.						LT	Recent	Neg	Invalid
15.						LT	Recent	Neg	Invalid
16.						LT	Recent	Neg	Invalid
17.						LT	Recent	Neg	Invalid
18.						LT	Recent	Neg	Invalid
19.						LT	Recent	Neg	Invalid
20.						LT	Recent	Neg	Invalid

Note: